

Crockmaster

COLOUR FASTNESS TO RUBBING TESTER 680

- > motorised
- > smooth, precision-engineered mechanism
- > pre-settable, digital counter
- > guaranteed accuracy of loading and stroke length
- > easy-to-use, polycarbonate specimen clamp
- > compliant with: EN ISO X12, AATCC 8 and M & S C8



Loading weight



Specimen clamp



Rubbing finger



Head arrangement



Digital counter

How to Order

902-415 CROCKMASTER – MODEL 680 *Without consumables and accessories*

CONSUMABLES BS/EN/ISO/M & S	ACCESSORIES	RECOMMENDED STARTING KITS
766-201 Grey Scale for assessing Staining 702-540 Crocking Cloths 50 x 50 mm – per pack (500) 701-217 Silicone Carbide Waterproof Paper 280 grit (M & S C8) – per pack (50)	794-998 Interchangeable Finger (Crockblock) for testing floor coverings and other pile fabrics (ISO 105 X12 and AATCC 165) 794-997 Interchangeable Token Holder, 50 Tokens and Weight for scuff testing of woven and knitted upholstery fabrics (BS 2543) 789-521 Tokens (BS 2543) – per pack (250)	BS/EN/ISO 902-415 1 x Crockmaster 766-201 1 x Grey Scale for assessing Staining 702-540 2 x packs (500) Crocking Cloths
CONSUMABLES AATCC 702-424 AATCC Style 3 Crocking Cloths 50 x 50 mm – per pack (1000) 766-510 AATCC Chromatic Transference Scale 9-step scale to evaluate colour transfer or staining – to facilitate evaluation by inexperienced personnel 766-513 AATCC Gray Scale for Staining	SPARE PARTS 521-690 Spare Clamp 117-499 Spare Counter 794-499 Spare Finger (16mm diameter) complete with Spring Clip 375-451 Spare Spring Clips – per pack (10)	M & S 902-415 1 x Crockmaster 766-201 1 x Grey Scale for assessing Staining 702-540 2 x Packs (500) Crocking Cloths 701-217 1 x Pack (50) Silicone Carbide Waterproof Paper
	CALIBRATION 201-255 ISO Certificate of Calibration for Crockmeter	AATCC 902-415 1 X Crockmaster 766-513 1 x AATCC Gray Scale for Staining 702-424 1 x Pack (1000) AATCC Style 3 Crocking Cloths



JAMES H. HEAL 
& CO. LTD. HALIFAX ENGLAND

Richmond Works, Halifax, West Yorkshire HX3 6EP, England.

Tel: +44 (0) 1422 366355. Fax: +44 (0) 1422 352440.

E-mail: info@james-heal.co.uk Internet: <http://www.james-heal.co.uk>

